

# FORUM

The Stieglmeyer-Group's company magazine



## Spotting TRENDS for the future with Burmeier

02. July 2017 // **Trade Fairs & Events**

Burmeier strives to the summit of a very successful year for trade fairs: At REHACARE in Düsseldorf, our home care specialists will be displaying products from October 4th to 7th that will set important trends for the upcoming years.

### **At REHACARE in Düsseldorf, exciting discoveries await**

**Burmeier strives to the summit of a very successful year for trade fairs: At REHACARE in Düsseldorf, our home care specialists will be displaying products from October 4th to 7th that will set important trends for the upcoming years.**

The machine-washable Dali-Wash care bed by Burmeier is best-suited for these tasks. At REHACARE, experience this bed together with the modern reprocessing system of RehaWash, our partner in cooperation. Our experts will tell you how easily and efficiently automated cleaning and disinfecting takes place.

The exceptional results will win you over.

Low-height beds for fall prevention and half-length safety sides for unhindered mobility also ensure high quality of life during home care. Our Regia care bed fulfils these criteria brilliantly while also being beautiful to look at and comfortable.

There is even more to discover at the Burmeier stand – the elegant, colourful Softcovers for our Dali beds or the comfortable bed inserts of the Lippe series, which can be fitted effortlessly into existing beds. Let us advise you in an inviting atmosphere.

**We will be happy to see you in Hall 5 at stand J24.**



Cosy and comfortable: the Westfalia IV care bed.

The Westfalia IV care bed, the newest model in the renowned Westfalia range, is a new addition to the Burmeier portfolio from this summer onwards, and is distinguished by its scissor-type frame which delivers great stability. Compared to the Westfalia III, its predecessor model, the design of the Westfalia IV has been advanced in many respects. The height adjustment range has been expanded downwards by 5 cm, giving it a range spanning 35 cm to 80 cm. This means occupants can always find a convenient height for getting in and out of bed, and nursing staff are able to work at a back-friendly height. The safe working load was increased to 185 kg, and the 24-volt drive system provides an optimised standard of safety. With a variety of mattress base sizes to choose from, the bed can be made to fit ideally into the available space, while the large selection of homelike head and footboards available allows the bed to blend in with its surroundings.